



Pure or not...?



WI SCAN

S y s t e m



Electronic Bowie & Dick *A new step in the WI SCAN System*

Innovation is our driving force

After the electronic sterilization controllers, Sterlab breaks new ground again by imagining a new electronic Bowie & Dick test.

A Universal Bowie & Dick test

The Sterlab Research and Development team has patented a new air and NCG detection concept to analyse autoclave steam, which brings a never yet obtained reliability in Bowie & Dick tests for all pre-vacuum autoclaves.

Integrated in the WI SCAN system

The WI SCAN electronic Bowie & Dick test will simplify controls, suppress the use of paper, optimise costs with the reader/software architecture of the existing WI SCAN System to read and save electronic data.

Conceived for long term development

As the WI SCAN electronic Bowie & Dick test can be reused year after year with no limit in time, the concept is far in advance of the single-use product, thus promoting responsible resource management and our contribution to the planet.

The WI SCAN electronic Bowie & Dick has been tested according to the EN ISO 11140-4 Norm.

Ref. WIBODI01 - Bowie & Dick test
Ref. WISCAN01 - WI SCAN Reader
Ref. WISOFT01 - Software
Ref. WICHAR01 - Charger for reader

