




SEALabel[®] et monoSEAL[®]

An unique concept...

The label holder folds and clips on the seal and is positioned on the front side of any container.

The oval window allows you to check if the **SEALabel[®]** is correctly locked.

The electronic traceability label is stuck on the **SEALabel[®]** holder.




Type 1 steam indicator in accordance with ISO 11140-1, turns from pink to black to check whether the container has already undergone a sterilization cycle or not.




OK turn from pink to black

References

monoSEAL[®] with ink indicator - Box of 1,000 units.



monoSEAL[®] without ink indicator - Box of 1,000 units.



SEALabel[®] + **monoSEAL[®]** with ink indicator.

Box of 1,000 units (500 **SEALabel[®]** + 500 **monoSEAL[®]**).



Patent pending

